
Optics Lab 3


Harsukh Singh

February 6, 2014

1 OVERVIEW OF EXPERIMENT

This lab involved the study and use of the Michelson Interferometer, the diagram of the interferometer is shown below. The Interferometer are usually used to study interference patterns of waves. Interferometers can usually be divided into two separate parts, the amplitude splitting and wave front splitting. The Michelson Interferometer is a amplitude splitting interfometer, so when a wave meets a boundary, a part of it is reflected and the other part is transmitted, both of these “new” waves have lower amplitudes compared to the original wave, so in a sense the amplitude has been split. If these new waves were brought together again it would create an interference pattern. The Michelson Interferometer does this by using two mirrors and a compensator, shown below the elements are drawn such that they appear in a straight line, for the sake of analysis.

Figure 1.1: Simplified diagram, with lightsource and mirror replaced by respective images formed by beam splitter.



With this diagram it is easy to see that the optical path difference is given by $OPD = 2d\cos(\theta)$. When the condition described by (1.1) is met, destructive interference will occur and an observer will see a circular fringe system.


$$2d\cos(\theta_m) = m\lambda_0 \quad (1.1)$$

In this equation d is the difference in path lengths of the two arms of the interferometer, m is an integer and θ_m is the angle of the m th dark ring. So the central dark fringe will be when $\theta_m = 0$ and can be represented by $2d = m_0\lambda_0$. To measure the wavelength of the source, one can move the arms of the mirror by a distance Δd and using the equation $\lambda = \frac{2\Delta d}{N}$.

2 COUNTING FRINGES

The set-up of the lab is shown below where, the laser is aligned to both mirrors on both arms of the interferometer so that the beams that meet the target are coinciding. This insures that beam is non-dispersive throughout the experiment and the interference pattern is observed. The translational stage is connected to a microcontroller for changing the displacement of the arms. Before the beam reaches the stage a -25 mm diverging lens into the beam to make things clearer and observe the fringe patterns better. Once this was done, translational stage

Figure 2.1: Experimental Setup



was moved and the fringes were counted, this is described by the table below.

$\lambda(nm)$	N	$\Delta d(m)$
658	48	0.0000158
649	49	0.0000159
615	53	0.0000163
644	50	0.0000161
620	50	0.0000155


(2.1)

The average was 637 nm which is close to the theoretical 633 nm.

3 VERIFYING THE TRANSLATIONAL STAGE VELOCITY

In this part of the lab, a photodetector was replaced as the target (Fig 2.1) and connected to an oscilloscope to measure the sinusoidal pattern of the fringes that appear on the central spot of the interfering waves, in an effort to verify the translational stage velocity. The result of sinusoidal fringe pattern is shown below,

Figure 3.1: Intensity measurement of Central Fringe with a translational stage velocity



From this figure and the change in distances we find that the stage velocity is $.07\text{mm/s}$, the new formula for the measurement of λ is,

$$\lambda = \frac{2v}{Nf} \quad (3.1)$$

where f is the frequency measured of the central fringe. With this measurement 640 nm.

4 PRECISE WAVELENGTH MEASUREMENT

In this part of the lab we verified two different situations, the first being the Fresnel- Arago Law, which says that two orthogonal, coherent P-states cannot interfere in the sense that no fringes result. This was the case when a polarizer was placed in front of the two mirrors of the arms of the interferometer and no fringes were found. The other one being the situation if a uncoated glass slide was placed into one of the arms, the result was a π phase shift because the beam has a lateral shift due to the glass but still emerges parallel to where it entered.